

## Human APEX1 protein (Recombinant) (Tag free) (STJP014463)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µm filtered solution in 10mM HEPES, 100mM KCl, 50% Glycerol, pH 7.4.
Molecular Weight	35.62 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Immunogen Region	Pro2-Leu318
------------------	-------------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--