

Human LEP/Leptin protein (Recombinant) (Tag free) (STJP014496)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM Tris-HCl, 50mM NaCl, 10% Trehalose, 0.02% Tween 80, pH 8.5.
Molecular Weight	16.1 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3952
Gene Symbol	LEP
UniProt ID	LEP_HUMAN
Immunogen Region	Val22-Cys167

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--