

Human AKR1C2 protein (Recombinant) (Tag free) (STJP014532)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 100mM NaCl, 1mM DTT, pH 8.0.
Molecular Weight	36.74 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	1646
Gene Symbol	AKR1C2
UniProt ID	AK1C2_HUMAN
Immunogen Region	Met1-Tyr323

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--