

Human SUMO3 protein (Recombinant) (N-His) (STJP014544)

GENERAL INFORMATION

| | |
|---------------|----------|
| Product Type | Proteins |
| Applications | ELISA/WB |
| Host / Source | E.coli |

PRODUCT PROPERTIES

| | |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation | Lyophilized from a 0.22 Mg/ml filtered solution in 20mM PB, 150mM NaCl, pH 7.2. |
| Molecular Weight | 13.8 kDa |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

TARGET INFORMATION

| | |
|------------------|-----------------------------|
| Gene ID | 6612 |
| Gene Symbol | SUMO3 |
| UniProt ID | SUMO3_HUMAN |
| Immunogen Region | Met1-Phe103 |

ADDITIONAL INFORMATION

| | |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|