

Human BID protein (Recombinant) (Tag free) (STJP014546)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM PB, 100mM KCl, pH 7.4.
Molecular Weight	21.99 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	637
Gene Symbol	BID
UniProt ID	BID_HUMAN
Immunogen Region	Met1-Asp195

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--