

Human LGALS4 protein (Recombinant) (C-His) (STJP014547)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mg/ml filtered solution in 50mM HEPES, 150mM NaCl, 5mM EDTA, 2mM TCEP, pH 7.4.
Molecular Weight	37.7 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3960
Gene Symbol	LGALS4
UniProt ID	LEG4_HUMAN
Immunogen Region	Met1-Ile323

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--