

Human CLDN18/Claudin-18 protein (Recombinant) (N-8*His) (STJP014548)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 25mM Tris-HCl, 25mM NaCl, 0.1% Triton X-100, 10% Glycerol, 0.1mM GSH, 0.01mM GSSG, pH 8.0.
Molecular Weight	19.7 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	Asp28-Leu76
------------------	-------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--