

Human UBE2H protein (Recombinant) (N-GST) (STJP014567)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.
Molecular Weight	47 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	7328
Gene Symbol	UBE2H
UniProt ID	UBE2H_HUMAN
Immunogen Region	Met1-Leu183

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--