

Human DYNLL1 protein (Recombinant) (N-His) (STJP014576)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 M filtered solution in 20mM Histidine-HCl, 10% Trehalose, 0.05% Tween 80, pH 6.0.
Molecular Weight	12.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	8655
Gene Symbol	DYNLL1
UniProt ID	DYL1_HUMAN
Immunogen Region	Met1-Gly89

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--