

Human CYCS protein (Recombinant) (C-His) (STJP014591)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM PB, 10% Trehalose, 200mM NaCl, 50% Glycerol, 0.05% Tween 80, pH 7.0.
Molecular Weight	12.8 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	54205
Gene Symbol	CYCS
UniProt ID	CYC_HUMAN
Immunogen Region	Gly2-Glu105

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--