

Human MC1R protein (Recombinant) (N-His) (STJP014596)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in PBS pH 7.4, 5% Trehalose, 2% Mannitol
Molecular Weight	14.42 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	4157
Gene Symbol	MC1R
UniProt ID	MSHR_HUMAN
Immunogen Region	[Met1-Glu37] & [Glu94-Asn118] & [Asp141-Arg163] & [Ala212-Ala240]

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--