

Human APBA1 protein (Recombinant) (N-His) (STJP014600)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli
PRODUCT PROPERTIES	
Formulation	Lyophilized from a solution in PBS pH 7.4, 5% Trehalose, 2% Mannitol.
Molecular Weight	21.22 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene ID	320
Gene Symbol	APBA1
UniProt ID	APBA1_HUMAN
Immunogen Region	Glu453-Pro623
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.