

## Human PRDX1 protein (Recombinant) (N-His & C-His) (STJP014604)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in PBS, 10% Glycerol, 0.1mM DTT, pH 6.0.
Molecular Weight	25.3 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">5052</a>
Gene Symbol	<a href="#">PRDX1</a>
UniProt ID	<a href="#">PRDX1_HUMAN</a>
Immunogen Region	Met1-Lys199

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--