

Human SMAD4 protein (Recombinant) (C-His) (STJP014621)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM Tris-HCl, pH 8.0.
Molecular Weight	61.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	4089
Gene Symbol	SMAD4
UniProt ID	SMAD4_HUMAN
Immunogen Region	Met1-Asp552

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--