

Human KPNB1 protein (Recombinant) (N-His) (STJP014632)

GENERAL INFORMATION		
	Product Type	Proteins
	Applications	ELISA/WB
	Host / Source	E.coli
PRODUCT PROPERTIES		
	Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 1mM DTT, 30% Glycerol, 100mM NaCl, pH 8.0.
	Molecular Weight	98.6 kDa
	Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION		
	Gene ID	3837
	Gene Symbol	KPNB1
	UniProt ID	IMB1_HUMAN
	Immunogen Region	Met1-Ala876 (Q559H)
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.