

## Human ACAT2 protein (Recombinant) (N-Trx-His) (STJP014686)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | ELISA/WB |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation         | Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 100mM NaCl, 1mM DTT, 10% Glycerol, pH 8.0.                                                               |
| Molecular Weight    | 59.1 kDa                                                                                                                                                              |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

### TARGET INFORMATION

|                  |                            |
|------------------|----------------------------|
| Gene ID          | <a href="#">39</a>         |
| Gene Symbol      | <a href="#">ACAT2</a>      |
| UniProt ID       | <a href="#">THIC_HUMAN</a> |
| Immunogen Region | Met1-Glu397                |

### ADDITIONAL INFORMATION

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|