

Human ITPA protein (Recombinant) (C-His) (STJP014687)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 5% Trehalose, 300mM NaCl, 30% Glycerol, 0.05% Tween 80, pH 8.0.
Molecular Weight	22.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3704
Gene Symbol	ITPA
UniProt ID	ITPA_HUMAN
Immunogen Region	Ala2-Ala194

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--