

Human LC3B/MAP1LC3B protein (Recombinant) (Tag free) (STJP014689)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM Citrate, 15% Trehalose, 0.05% Tween 80, 2mM DTT, pH 4.0.
Molecular Weight	14.8 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	81631
Gene Symbol	MAP1LC3B
UniProt ID	MLP3B_HUMAN
Immunogen Region	Met1-Val125

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--