

## Human ATG5 protein (Recombinant) (Tag free) (STJP014693)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 0.2M NaCl, 1mM DTT, 40% Glycerol, pH 8.0.
Molecular Weight	32.4 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">9474</a>
Gene Symbol	<a href="#">ATG5</a>
UniProt ID	<a href="#">ATG5_HUMAN</a>
Immunogen Region	Met1-Asp275

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	----------------------------------------------------------------------------------------------------------------------