

Human CCL26/Eotaxin-3/TSC-1 protein (Recombinant) (Tag free) (STJP014708)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM PB, 150mM NaCl, pH 7.2.
Molecular Weight	8.21 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	10344
Gene Symbol	CCL26
UniProt ID	CCL26_HUMAN
Immunogen Region	Ser27-Leu94

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--