

## Human CCL26/Eotaxin-3/TSC-1 protein (Recombinant) (Tag free) (STJP014709)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 1mM EDTA, 20% Glycerol, pH 9.0.
Molecular Weight	8.53 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">10344</a>
Gene Symbol	<a href="#">CCL26</a>
UniProt ID	<a href="#">CCL26_HUMAN</a>
Immunogen Region	Thr24-Leu94

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--