

Human CCL27/CTACK protein (Recombinant) (Tag free) (STJP014710)

GENERAL INFORMATION

| | |
|---------------|----------|
| Product Type | Proteins |
| Applications | ELISA/WB |
| Host / Source | E.coli |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Formulation | Lyophilized from a 0.22 µm filtered solution in 20mM PB, 100mM NaCl, 6% Trehalose, 4% Mannitol, 0.05% Tween 80, pH 7.0. |
| Molecular Weight | 10.1 kDa |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

TARGET INFORMATION

| | |
|------------------|-----------------------------|
| Gene ID | 10850 |
| Gene Symbol | CCL27 |
| UniProt ID | CCL27_HUMAN |
| Immunogen Region | Phe25-Gly112 |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|