

Human DPP3 protein (Recombinant) (N-His) (STJP014733)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 25mM Tris-HCl, 10mM NaCl, 2mM DTT, 0.05% Tween 80, 20% Glycerol, pH 8.0.
Molecular Weight	84.8 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	10072
Gene Symbol	DPP3
UniProt ID	DPP3_HUMAN
Immunogen Region	Met1-Ala737

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--