

Human GPER1 protein (Recombinant) (N-His-KSI) (STJP014746)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli
PRODUCT PROPERTIES	
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	27.65 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene ID	2852
Gene Symbol	GPER1
UniProt ID	GPER1_HUMAN
Immunogen Region	[Asp2-Ser62] & [Gly328-Val375]
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.