

Human CCDC111/PRIMPOL protein (Recombinant) (C-His) (STJP014771)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% SKL, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	25.36 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	201973
Gene Symbol	PRIMPOL
UniProt ID	PRIPO_HUMAN
Immunogen Region	Arg3-Sea202

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--