

Human ACY3 protein (Recombinant) (N-His) (STJP014801)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM PB, 10% Trehalose, 10% Glycerol, 1 mM DTT, 0.02% Tween 80, pH7.6.
Molecular Weight	37.4 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	91703
Gene Symbol	ACY3
UniProt ID	ACY3_HUMAN
Immunogen Region	Met1-Ser319

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--