

Human NECAP2 protein (Recombinant) (N-His & C-His) (STJP014819)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in PBS, 2mM EDTA, 1mM DTT, pH 7.4.
Molecular Weight	31.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	55707
Gene Symbol	NECAP2
UniProt ID	NECP2_HUMAN
Immunogen Region	Met1-Phe263

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--