

Human NIP7 protein (Recombinant) (N-His) (STJP014825)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mg/ml filtered solution in 20mM Tris-HCl, 100mM NaCl, pH 8.0.
Molecular Weight	22.6 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	51388
Gene Symbol	NIP7
UniProt ID	NIP7_HUMAN
Immunogen Region	Met1-Thr180

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--