

Human PPIL1 protein (Recombinant) (N-His) (STJP014828)

GENERAL INFORMATION

| | |
|---------------|----------|
| Product Type | Proteins |
| Applications | ELISA/WB |
| Host / Source | E.coli |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Formulation | Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, pH 8.0. |
| Molecular Weight | 20.4 kDa |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

TARGET INFORMATION

| | |
|------------------|-----------------------------|
| Gene ID | 51645 |
| Gene Symbol | PPIL1 |
| UniProt ID | PPIL1_HUMAN |
| Immunogen Region | Met1-Gly166 |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|