

## Human SMDT1 protein (Recombinant) (N-GST & C-His) (STJP014831)

| GENERAL INFORMATION    |                                                                                                                                                                           |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type           | Proteins                                                                                                                                                                  |
| Applications           | ELISA/WB                                                                                                                                                                  |
| Host / Source          | E.coli                                                                                                                                                                    |
| PRODUCT PROPERTIES     |                                                                                                                                                                           |
| Formulation            | Lyophilized from a solution in PBS pH 7.4, 0.02% SKL, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                |
| Molecular Weight       | 39.44 kDa                                                                                                                                                                 |
| Storage Instruction    | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |
| TARGET INFORMATION     |                                                                                                                                                                           |
| Gene ID                | <a href="#">91689</a>                                                                                                                                                     |
| Gene Symbol            | <a href="#">SMDT1</a>                                                                                                                                                     |
| UniProt ID             | <a href="#">EMRE_HUMAN</a>                                                                                                                                                |
| Immunogen Region       | Met1-Asp107                                                                                                                                                               |
| ADDITIONAL INFORMATION |                                                                                                                                                                           |
| Note                   | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.                                                      |