

CHIKV Spike glycoprotein E2 protein (Recombinant) (N-GST) (STJP014932)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli
PRODUCT PROPERTIES	
Formulation	Lyophilized from a solution in PBS pH 7.4, 5% Trehalose, 2% Mannitol.
Molecular Weight	65.06 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Immunogen Region	Asn330-Gln666
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.