

Dog CD68/Macrosialin protein (Recombinant) (C-Fc) (STJP015098)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/WB/ELISA/Immunogen
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	57.13 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -147°C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	Ser22-Asn295
------------------	--------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--