

Mouse CD69 protein (Recombinant) (C-Fc) (STJP015130)

GENERAL INFORMATION	
Product Type	Proteins
Applications	SDS-PAGE/WB/ELISA/Immunogen
Host / Source	Mammalian Cells
PRODUCT PROPERTIES	
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	44.16 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -179°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene ID	12515
Gene Symbol	Cd69
UniProt ID	CD69_MOUSE
Immunogen Region	Asn62-Arg199
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.