

## Human APOBEC3C protein (Recombinant) (N-His) (STJP015165)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Lyophilized from PBS pH7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	23.71 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -214°C for twelve months from the date of receipt.

### TARGET INFORMATION

Immunogen Region	Met12-Gln190
------------------	--------------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--