

Human CCNB1 protein (Recombinant) (N-GST & C-His) (STJP015178)

GENERAL INFORMATION

| | |
|---------------|-----------------------------|
| Product Type | Proteins |
| Applications | ELISA/Immunogen/SDS-PAGE/WB |
| Host / Source | E.coli |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Formulation | Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol. |
| Molecular Weight | 58.44 kDa |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -227°C for twelve months from the date of receipt. |

TARGET INFORMATION

| | |
|------------------|-----------------------------|
| Gene ID | 891 |
| Gene Symbol | CCNB1 |
| UniProt ID | CCNB1_HUMAN |
| Immunogen Region | Val171-Val433 |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|