

Human APOBEC3A protein (Recombinant) (N-His) (STJP015194)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from PBS pH7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	25.17 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -243°C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	Met1-Asn199
------------------	-------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--