

Human ERCC8 protein (Recombinant) (N-His) (STJP015905)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	43.36 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	1161
Gene Symbol	ERCC8
UniProt ID	ERCC8_HUMAN
Immunogen Region	Met1-Pro370

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	----------------------------------------------------------------------------------------------------------------------