

## Human SLC25A1 protein (Recombinant) (N-GST & C-His) (STJP016092)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	44.78 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">6576</a>
Gene Symbol	<a href="#">SLC25A1</a>
UniProt ID	<a href="#">TXTP_HUMAN</a>
Immunogen Region	Glu47-Arg85&Glu106-Gly128&Glu144-Gln182&Leu241-Gly278

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--