

Stem cell factor receptor protein (STJP016396)

GENERAL INFORMATION

Product Type	Proteins
Applications	Block
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	When reconstituted in 0.5 mL sterile phosphate-buffered saline, the solution will contain 1% human serum albumin (HSA) and 10% trehalose.
Purification	Purified. The SCFR-Fc Chimera consists of 5-45% carbohydrate by weight.
Storage Instruction	It is recommended that 0.5 mL of sterile phosphate-buffered saline be added to the vial. Lyophilized products should be stored at 2 to 8°C. Following reconstitution short-term storage at 2-8°C is recommended, and longer-term storage of aliquots at -18

TARGET INFORMATION

Gene ID	3815
Gene Symbol	KIT
UniProt ID	KIT_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---