

L-selectin protein (STJP016398)

GENERAL INFORMATION

Product Type	Proteins
Applications	Block
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	When reconstituted in 0.5 mL sterile phosphate-buffered saline, the solution will contain 1% human serum albumin (HSA) and 10% trehalose.
Purification	Purified. The L-Selectin-Fc chimera consists of 20-45% carbohydrate by weight.
Storage Instruction	Reconstitute with sterile-filtered phosphate-buffered saline to a concentration of 100 µg/mL. Lyophilized products should be stored at 2 to 8°C. Following reconstitution, short-term storage at 2-8°C is recommended and longer-term storage of aliquots

TARGET INFORMATION

Gene ID	6402
Gene Symbol	SELL
UniProt ID	LYAM1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---