

## Casein, beta A2 protein (STJP016402)

### GENERAL INFORMATION

Product Type	Proteins
Applications	WB
Host / Source	Natural A2/A2 Cows

### PRODUCT PROPERTIES

Formulation	Protein is provided as 50 $\mu$ L aliquot containing protease inhibitor at 0.4 mg/mL in TBS, pH 8.5. Does not contain preservatives. Use sterile methods when handling. Tris-based buffers are subject to pH changes with temperature. Keep protein solutions a
Purification	Purified
Storage Instruction	Protein is shipped on ice packs. Upon receipt, divide into aliquots and store product long-term at -10°C to -20°C. Do not store at -80°C. Purified proteins will precipitate from the solution when cold or at a pH below 8.0. Purified Caseins are very sens

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---