

Parkinson disease protein 2 protein (STJP016404)

GENERAL INFORMATION

Product Type	Proteins
Applications	Block
Host / Source	Synthetic

PRODUCT PROPERTIES

Purification	Purified
Storage Instruction	This peptide is soluble in PBS and water. After reconstitution keep at -20°C. Avoid repetitive freeze/thaw cycles.

TARGET INFORMATION

Gene ID	5071
Gene Symbol	PRKN
UniProt ID	PRKN_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---