

Human Adiponectin/Acrp30 protein (Recombinant) (No-Tag) (STJP016408)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 1mM DTT
Molecular Weight	25.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9370
Gene Symbol	ADIPOQ
UniProt ID	ADIPO_HUMAN
Immunogen Region	15-244aa
Immunogen Sequence	MGHDQETTTQ GPGVLLPLPK GACTGWMAGI PGHPGHNGAP GRDGRDGTGP EKGEKGDPLG IGPKGDIGET GVPGAEGPRG FPGIQGRKGE PGEGAYVYRS AFSVGLETYV TIPNMPIRFT KIFYNQNNHY DGSTGKFHCN IPGLYYFAYH ITVYMKDVKV SLFKKDKAML FTYDQYQENN VDQASGSVLL HLEVGQVWL QVYGEGERNG LYADNDNDST FTGFLLYH

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human Adiponectin Acrp30 protein 3µg under reducing conditions with Coomassie blue staining