

Human AK4 (C22S) protein (Recombinant) (His-Tag) (STJP016411)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.5) containing 2mM DTT, 30% Glycerol, 0.1M NaCl
Molecular Weight	29.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	205
Gene Symbol	AK4
UniProt ID	KAD4_HUMAN
Immunogen Region	1-223aa
Immunogen Sequence	MASK LLRAVILGPP GSGKGTVSQR IAQNFGLQHL SSGHFLRENI KASTEVGEMA KQYIEKSLLV PDHVITRLMM SELENRRGQH WLLDGFPRTL GQAEALDKIC EVDLIVSLNI PFETLKDRLS RRWIHPPSGR VYNLDFNPPH VHGIDDVTGE PLVQQEDDKP EAVAARLRQY KDVAKPVIEL YKSRGV

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human AK4 C22S protein under reducing conditions stained with Coomassie blue