

## Human AHA1 protein (Recombinant) (No-Tag) (STJP016414)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 2mM EDTA
Molecular Weight	36.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">10598</a>
Gene Symbol	<a href="#">AHSA1</a>
UniProt ID	<a href="#">AHSA1_HUMAN</a>
Immunogen Region	19-337aa
Immunogen Sequence	MATNVNNWHW TERDASNWST DKLKTFLAV QVQNEEGKCE VTEVSKLDGE ASINNRKGKL IFFYEWSVKL NWTGTSKSGV QYKGHVEIPN LSDENSVDV EISVSLAKDE PDTNLVALMK EEGVKLLREA MGIYISTLKT EFTQGMILPT MNGESVDPVG QPALKTEERK AKPAPSKTQA RPPGVKIPTC KITLKETFLT SPEELYRVFT TQELVQAFTH APATLEAD

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human AHA1 protein under reducing conditions Coomassie blue stain