

Human Latexin protein (Recombinant) (No-Tag) (STJP016426)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 50mM NaCl, 10% Glycerol
Molecular Weight	25.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56925
Gene Symbol	LXN
UniProt ID	LXN_HUMAN
Immunogen Region	1-222aa
Immunogen Sequence	MEIPPTNYP SRAALVAQNY INYQQGTPHR VFEVQKVKQA SMEDIPGRGH KYRLKFAVEE IIQKQVKVNC TAEVLYPSTG QETAPEVNFET FEGETGKNPD EEDNTFYQRL KSMKEPLEAQ NIPDNFGNVS PEMTLVLHLA WVACGYIIWQ NSTEDTWYKM VKIQTVMQVQ RNDDFIELDY TILLHNIASQ EIIPWQMQL WHPQYGTKVK HNSRLPKEVQ LE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human Latexin protein 3 µg under reducing conditions stained with Coomassie blue