

Human Annexin A1/ANXA1 protein (Recombinant) (No-Tag) (STJP016440)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 100mM NaCl, 1mM DTT, 10% Glycerol
Molecular Weight	38.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	301
Gene Symbol	ANXA1
UniProt ID	ANXA1_HUMAN
Immunogen Region	1-346aa
Immunogen Sequence	MAMVSEFLKQ AWFIEEEQE YVQTVKSSKG GPGSAVSPYP TFPSSDVAA LHKAIMVKGV DEATIIDLTKRNNAQRQQI KAAYLQETGK PLDETLKKAL TGHLEEVVLA LLKTPAQFDA DELRAAMKGL GTDEDTLIEI LASRTNKEIR DINRVYREEL KRDLAKDITS DTSGDFRNAL LSLAKGDRSE DFGVNEDLAD SDARALYEAG ERRKGTDVNV FNTILTTR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant human Annexin A1 ANXA1 protein
SDS-PAGE gel 3µg reducing conditions
Coomassie blue stain