

Human FABP3/H-FABP protein (Recombinant) (His-Tag) (STJP016441)

GENERAL INFORMATION

| | |
|---------------|----------|
| Product Type | Proteins |
| Applications | SDS-PAGE |
| Host / Source | E.coli |

PRODUCT PROPERTIES


| | |
|---------------------|---|
| Concentration | 1 mg/mL |
| Formulation | Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol |
| Molecular Weight | 19.1 kDa |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|--------------------|---|
| Gene ID | 2170 |
| Gene Symbol | FABP3 |
| UniProt ID | FABPH_HUMAN |
| Immunogen Region | 1-133aa |
| Immunogen Sequence | MVD AFLGTWKLVD SKNFDDYMKV LGVGFATRQV ASMTKPTTII EKNGDILTLK THSTFKNTEI SFKLGVEFDE TTADDRKVKS IVTLDGGKLV HLQKWGGQET TLVRELIDGK LILTLTHGTA VCTRTRYEKEA |

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human FABP3 H-FABP protein under reducing conditions stained with Coomassie blue