

Human GRO beta/CXCL2 protein (Recombinant) (His-Tag) (STJP016446)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0)
Molecular Weight	10.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2920
Gene Symbol	CXCL2
UniProt ID	CXCL2_HUMAN
Immunogen Region	35-107aa
Immunogen Sequence	APLATELRRC QCLQLQGIH LKNIQSVKVK SPGPHCAQTE VIATLKNGQK ACLNPASPMV KKIIEKMLKN GKSN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human GRO beta CXCL2 protein under reducing conditions with Coomassie blue stain