

Human FKBP52/FKBP4 protein (Recombinant) (His-Tag) (STJP016450)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Enzyme Activity
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	53.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2288
Gene Symbol	FKBP4
UniProt ID	FKBP4_HUMAN
Immunogen Region	1-459aa
Immunogen Sequence	MTAEEMKATE SGAQSAPLPM EGVDISPKQD EGVLKVIKRE GTGTEMPMIG DRVHVHYTGW LLDGTFKFDSS LDRKDKFSFD LGKGVEIKAW DIAIATMKVG EVCHITCKPE YAYGSAGSPP KIPPNATLVF EVELFEFKGE DLTEEDGGI IRRIQTRGEG YAKPNEGAIV EVALEGYYKD KLFDQRELR F EIGEGENLDL PYGLER

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human FKBP52/FKBP4 protein 3µg under reducing conditions
Coomassie blue stain